

A mobile workstation consisting of a white base unit on four wheels, a white table with two tablets, and a yellow control unit on top. Red cables connect the control unit to the tablets. The workstation is in a classroom with a green floor and white walls. Another similar workstation is visible in the background.

Hi, Ich bin der
medienBuddy®!



ALLE VORTEILE AUF EINEN BLICK

Intuitiv

Begeisterendes Design, einfach zu bedienen, weckt die Freude an Technik bei allen Schülern

Sicher

Gas- und Stromversorgung durch Freigabe der Lehrer mit RFID-Chip, automatische Zentralverriegelung der Rollen, Not-Aus innerhalb direkter Reichweite, Neigungs- und Erschütterungssensor

Ergonomisch

Steckdosen und Gasanschlüsse auf Tischhöhe, keine langen Kabelverbindungen nötig, Medienversorgung direkt am Tisch

Wissenschaftlich

Keine aufwändige Projektplanung nötig, keine Deckenkonstruktionen oder Bodenarbeiten erforderlich, einfache Stromversorgung über Standard-Steckdose

Flexibel

Aufstellen, anschliessen und loslegen – drinnen oder draussen!

Nachhaltig

Langlebiger Lithium-Eisen-Phosphat-Akku, frei von Kobalt und Blei, zu 100 % recyclebar