


Hello, I'm
medienBuddy®!





info@waldner-benelux.nl

THE ADVANTAGES AT A GLANCE

Intuitive

Emotional appeal, engaging product design, easy to operate, promotes enthusiasm for technology among all students

Safe

Gas and power supply with teacher's control via RFID chip, automatic central locking of the castors, emergency shutdown within direct access range, tilt and shock sensor

Ergonomic

Sockets and gas connections at table height, no long cable connections necessary, media supply directly at the table

Economical

No elaborate project planning required, no structural ceiling or floor work necessary on site, easy to connect – a standard power socket is sufficient

Flexible

Set up, connect and you're ready to go – inside or outside!

Sustainable

Long-lasting lithium iron phosphate battery, cobalt- and lead-free, 100% recyclable